

LAMP2 / CD107b Mouse anti-Human Monoclonal (FITC) (H4B4) Antibody - LS-C86765 - LSBio	
<b>CatalogID:</b>	LS-C86765
<b>Target:</b>	lysosomal-associated membrane protein 2 (LAMP2)
<b>Synonyms:</b>	LAMP2 Antibody, CD107b Antibody, CD107b antigen Antibody, LAMPB Antibody, LGP110 Antibody, LAMP-2 Antibody
<b>Family / Subfamily:</b>	Cancer
<b>Host</b>	LAMP2 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	H4B4
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	LAMP2 / CD107b antibody was raised against Human
<b>Reactivity:</b>	Human
<b>Presentation:</b>	PBS, sodium azide
<b>Recommended Storage:</b>	Store at 4°C. Do not freeze.
<b>Uses:</b>	Immunofluorescence (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	