

PECAM-1 / CD31 Rat anti-Mouse Monoclonal (FITC) (390) Antibody - LS-C86798 - LSBio

<b>CatalogID:</b>	LS-C86798
<b>Target:</b>	platelet/endothelial cell adhesion molecule 1 (PECAM1)
<b>Synonyms:</b>	PECAM1 Antibody, CD31 Antibody, CD31 antigen Antibody, EndoCAM Antibody, PECA1 Antibody, PECAM-1 Antibody, GPIIA Antibody
<b>Host</b>	PECAM1 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Clone Name:</b>	390
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	PECAM-1 / CD31 antibody was raised against Mouse
<b>Specificity:</b>	This antibody is specific for Mouse CD31/PECAM-1, Mr 130-140 kD.
<b>Reactivity:</b>	Mouse
<b>Presentation:</b>	PBS, sodium azide
<b>Recommended Storage:</b>	Store at 4°C. Do not freeze.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry, Blocking (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences