

CD8A / CD8 Alpha Rat anti-Mouse Monoclonal (PE) (53-6.7) Antibody - LS-C86825 - LSBio	
CatalogID:	LS-C86825
Target:	CD8a molecule
Synonyms:	CD8A Antibody, CD8a molecule Antibody, CD8a antigen Antibody, OKT8 T-cell antigen Antibody, Leu2 T-lymphocyte antigen Antibody, T8 T-cell antigen Antibody, T cell co-receptor Antibody, p32 Antibody, CD8 Antibody, Leu2 Antibody, T-cell antigen Leu2 Antibody
Host	CD8A antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	53-6.7
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD8A / CD8 Alpha antibody was raised against Mouse
Immunogen:	CD8A / CD8 Alpha antibody was raised against mouse thymus or spleen.
Reactivity:	Mouse
Presentation:	PBS, sodium azide
Recommended Storage:	Store at 4°C. Product is photosensitive and should be protected from light.
Uses:	IHC, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	