

THY1 / CD90 Mouse anti-Human Monoclonal (PE) (F15-42-1) Antibody - LS-C86828 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C86828 |
| Target: | Thy-1 cell surface antigen (THY1) |
| Synonyms: | THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody |
| Host | THY1 antibody was produced in Mouse |
| Clonality: | Monoclonal |
| Isotype: | IgG1 |
| Clone Name: | F15-42-1 |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | THY1 / CD90 antibody was raised against Human |
| Reactivity: | Human |
| Purification: | Tissue culture supernatant |
| Presentation: | Culture supernatant diluted in PBS, 10 mM sodium azide, 1 mg/ml BSA |
| Recommended Storage: | Store at 4°C. Product is photosensitive and should be protected from light. |
| Uses: | Undetermined |
| Size: | 100 tst |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences