

Vimentin Mouse anti-Bovine Monoclonal (FITC) (VIM 3B4) Antibody - LS-C86895 - LSBio

CatalogID:	LS-C86895
Target:	vimentin
Synonyms:	VIM Antibody, Vimentin Antibody
Host	VIM antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	VIM 3B4
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Vimentin antibody was raised against Bovine
Antigen Type:	Purified protein
Immunogen:	Vimentin antibody was raised against vimentin purified from bovine lens.
Specificity:	Specific for the intermediate filament protein vimentin, tested in Human, Bovine, Dog, Chicken, Monkey and amphibians
Reactivity:	Bovine
Purification:	Protein A purified
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
Uses:	Undetermined
Size:	250 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences