

MYC / c-Myc Chicken anti-Human Polyclonal Antibody - LS-C86899 - LSBio

| | |
|-----------------------------|--|
| CatalogID: | LS-C86899 |
| Target: | v-myc avian myelocytomatosis viral oncogene homolog |
| Synonyms: | MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody |
| Host | MYC antibody was produced in Chicken |
| Clonality: | Polyclonal |
| Immunogen Species: | MYC / c-Myc antibody was raised against Human |
| Antigen Type: | Synthetic peptide - KLH conjugated |
| Immunogen: | MYC / c-Myc antibody was raised against synthetic peptide (hu c-myc epitope - KLH). |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.1% sodium azide |
| Recommended Storage: | Store at 4°C. Do not freeze. |
| Uses: | Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences