

CD45R Mouse anti-Human Monoclonal Antibody - LS-C88951 - LSBio	
<b>CatalogID:</b>	LS-C88951
<b>Target:</b>	CD45R
<b>Host</b>	CD45R antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Immunogen Species:</b>	CD45R antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD45R antibody was raised against human lymphocytes.
<b>Specificity:</b>	Recognizes human CD45 (LCA).
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Concentration:</b>	1 mg/ml
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	