

MCAM / CD146 Mouse anti-Chicken Monoclonal (PE) Antibody - LS-C8933 - LSBio

CatalogID:	LS-C8933
Target:	melanoma cell adhesion molecule (MCAM)
Synonyms:	MCAM Antibody, CD146 Antibody, CD146 antigen Antibody, Gicerin Antibody, Melanoma adhesion molecule Antibody, MUC18 Antibody, p1h12 Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	MCAM antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MCAM / CD146 antibody was raised against Chicken
Immunogen:	MCAM / CD146 antibody was raised against pooled Ed15 and Ed16 thymocytes from H. B15, H. B19 and H. B21 chicken embryos.
Specificity:	Recognizes chicken MCAM, Mr 98kD (protein core ~84kD).
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. Immunohistochemistry (acetone-fixed, frozen tissue sections). Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences