

AP50 / AP2M1 Rabbit anti-Human Polyclonal (aa150-200) Antibody - LS-C92140 - LSBio

CatalogID:	LS-C92140
Target:	adaptor-related protein complex 2, mu 1 subunit (AP2M1)
Synonyms:	AP2M1 Antibody, Adaptin-mu2 Antibody, AP-2 complex subunit mu Antibody, AP-2 mu 2 chain Antibody, AP-2 mu chain Antibody, AP50 Antibody, CLAPM1 Antibody, Mu2 Antibody, KIAA0109 Antibody, HA2 50 kDa subunit Antibody
Host	AP2M1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	AP50 / AP2M1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP50 / AP2M1 antibody was raised against synthetic peptide. within residues 150-200 of human AP2M1 protein (KLH) (Swiss Prot: Q96CW1).
Specificity:	Recognizes human AP2M1 protein.
Epitope:	aa150-200
Reactivity:	Human
Purification:	Immunoaffinity purified
Reconstitution:	Sterile 40-50% glycerol
Presentation:	PBS, pH 7.4, 0.02% sodium azide
Recommended Storage:	Store lyophilized powder -20°C. Reconstitute with sterile 40-50% glycerol, aliquot and store at -20°C.
Usage Summary:	Suitable for use in ELISA and Western Blot. ELISA: 0.05-0.2 ug/ml. Western Blot: 0.5-1 ug/ml.
Uses:	Western blot (0.5 - 1 µg/ml), ELISA (0.05 - 0.2 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	40 µg
Concentration:	0.5 mg/ml (after resuspension)
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences