

BrdU Mouse Monoclonal Antibody - LS-C92198 - LSBio

CatalogID:	LS-C92198
Target:	Bromodeoxyuridine (BrdU)
Host	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Specificity:	Recognizes Bromodeoxyuridine.
Purification:	Protein A purified
Presentation:	10 mM PBS, pH 7.4, 150 mM sodium chloride, 0.1% sodium azide
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry and Immunohistochemistry.
Uses:	IHC, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences