

GGA3 Rabbit anti-Human Polyclonal Antibody - LS-C94190 - LSBio

CatalogID:	LS-C94190
Target:	golgi-associated, gamma adaptin ear containing, ARF binding protein 3 (GGA3)
Synonyms:	GGA3 Antibody, KIAA0154 Antibody
Host	GGA3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	GGA3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GGA3 antibody was raised against a synthetic peptide from the VHS domain of human GGA3 (Golgi-localized, gamma ear-containing, ARF-binding protein 3) conjugated to an immunogenic carrier protein was used as the immunogen. The antigen is homologous in rat, mouse, chicken.
Specificity:	Appears to be specific for GGA3.
Reactivity:	Human, Mouse, Rat
Purification:	Purified
Reconstitution:	Reconstitute with sterile water. Possible additional volumes for resuspension: 500 μ l
Presentation:	Lyophilized. Centrifuge to remove any insoluble material
Recommended Storage:	Maintain the lyophilized/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Usage Summary:	IHC: Antigen retrieval is essential for use on paraffin sections.
Uses:	Western blot (10 - 50 μ g/ml) (Optimal dilution to be determined by the researcher)
Size:	500 μ g
Concentration:	1 mg/ml (after resuspension)
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences