

Complement C1QA Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C95297 - LSBio

CatalogID:	LS-C95297
Target:	complement component 1, q subcomponent, A chain (C1QA)
Synonyms:	C1QA Antibody, Complement C1QA Antibody
Host	C1QA antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Complement C1QA antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	Complement C1QA antibody was raised against c1q isolated from the normal human serum.
Specificity:	Reacts with human C1q complement. Traces of contaminating antibodies have been removed by solid phase absorption with human plasma proteins. This antibody is well suited for tissue section staining but may be used for other applications also.
Reactivity:	Human
Purification:	Ig fraction
Presentation:	0.05% sodium azide
Recommended Storage:	Store at 2°C to 8°C. Do not freeze. Because product is photosensitive, store in a dark location.
Uses:	Undetermined
Size:	1 ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences