

CDH1 / E Cadherin Mouse Monoclonal (67A4) Antibody - LS-C96302 - LSBio	
CatalogID:	LS-C96302
Target:	cadherin 1, type 1, E-cadherin (epithelial) (CDH1)
Synonyms:	CDH1 Antibody, CD324 Antibody, CDHE Antibody, Cell-CAM 120/80 Antibody, E Cadherin Antibody, E-Cadherin Antibody, ECAD Antibody, Epithelial cadherin Antibody, Arc-1 Antibody, UVO Antibody, Cadherin-1 Antibody, CAM 120/80 Antibody, CD324 antigen Antibody, LCAM Antibody, Uvomorulin Antibody
Family / Subfamily:	Cadherin / not assigned-Cadherin
Host	CDH1 antibody was produced in Mouse
Clonality:	Monoclonal
Clone Name:	67A4
Immunogen:	CDH1 / E Cadherin antibody was raised against breast carcinoma cell line T47D
Specificity:	Recognizes E-Cadherin, also called L-CAM (Liver Cell Adhesion Molecule).
Reactivity:	Human
Purification:	Affinity purified
Reconstitution:	0.5ml distilled water
Presentation:	Lyophilized from PBS pH 7.2, 10mg/ml BSA, 0.1% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	IHC - Frozen (1:50 - 1:100) (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 10/2/2014 © 2014 LifeSpan BioSciences	