

THY1 / CD90 Mouse Monoclonal (OX-7) Antibody - LS-C96412 - LSBio	
<b>CatalogID:</b>	LS-C96412
<b>Target:</b>	Thy-1 cell surface antigen (THY1)
<b>Synonyms:</b>	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
<b>Host</b>	THY1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Clone Name:</b>	OX-7
<b>Specificity:</b>	Recognizes the 18kDa Thy 1.1 antigenic determinant on rat and mouse cells.
<b>Reactivity:</b>	Rat
<b>Purification:</b>	Purified
<b>Reconstitution:</b>	0.5ml distilled water
<b>Presentation:</b>	Lyophilized from PBS pH 7.2, 0.02% sodium azide
<b>Recommended Storage:</b>	Store at 4°C.
<b>Uses:</b>	IHC - Frozen (1:200) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	