



POU5F1 / OCT4 Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-C97710 - LSBio	
CatalogID:	LS-C97710
Target:	POU class 5 homeobox 1 (POU5F1)
Synonyms:	POU5F1 Antibody, Oct-4 Antibody, OCT3 Antibody, OCT4 Antibody, OTF4 Antibody, POU class 5 homeobox 1 Antibody, Octamer-binding protein 3 Antibody, Octamer-binding protein 4 Antibody, OTF-3 Antibody, Oct-3 Antibody, OTF3 Antibody
Host	POU5F1 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	POU5F1 / OCT4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	POU5F1 / OCT4 antibody was raised against kLH conjugated synthetic peptide selected from the N-terminal region of human 39724.
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Ammonium sulfate precipitation
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	Immunofluorescence (1:50 - 1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Uses Not Recommended:	IHC - Paraffin
Size:	400 µl

Immunofluorescence Image:



Immunofluorescence of OCT3(OCT4) antibody (N-term) LS-C97710 in HeLa cells. 0.025 mg/ml primary antibody was followed by Alexa-Fluor-546-conjugated donkey anti-rabbit IgG (H+L). Alexa-Fluor-546 emits orange fluorescence. Blue counterstaining is DAPI.

Western Blot Image:



The anti-Oct3 N-term Pab LS-C97710 is used in Western blot to detect Oct3 in A375 cell lysate.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/3/2014

© 2014 LifeSpan BioSciences