

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

Product-No.:	01-19-203
Product-Name:	micromer ®
surface:	streptavidin
size:	2 µm
solid content:	25 mg/ml
quantity:	1 ml
Standard deviation:	< 5 % (C.V.)
Shape:	spherical
Density:	1.03 g/ccm
Stable in:	aqueous buffers
Not stable in:	organic solvents
Product form:	suspension in PBS buffer (pH = 7.4), stabilized in 0.02 % sodium azide
Particles per ml:	6.0*10E9
Particles per mg:	2.3*10E8
Protein binding capacity:	900 ng streptavidin / mg particle
Additional remarks:	Storage at 2 - 8°C for 3 months