

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

Product-No.:	01-20-102
Product-Name:	micromer ®
surface:	proteinA
size:	100 nm
solid content:	10 mg/ml
quantity:	1 ml
Polydispersity index:	< 0.2
Shape:	spherical
Density:	1.03 g/ccm
Stable in:	aqueous buffers, methanol, ethanol, DMSO
Not stable in:	organic solvents
Product form:	suspension in PBS buffer (pH = 7.4), stabilized in 0.02 % sodium azide
Particles per ml:	1.9*10E13
Particles per mg:	1.9*10E12
Additional remarks:	Storage at 2 - 8°C for 3 months