

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

| | |
|---------------------|---------------------------------|
| Product-No.: | 01-25-303 |
| Product-Name: | micromer ® |
| surface: | collagen |
| size: | 3 µm |
| solid content: | 25 mg/ml |
| quantity: | 5 ml |
| Standard deviation: | < 5 % (C.V.) |
| Shape: | spherical |
| Density: | 1.03 g/ccm |
| Stable in: | aqueous buffers |
| Not stable in: | organic solvents |
| Product form: | aqueous suspension |
| Particles per ml: | 1.7*10E9 |
| Particles per mg: | 6.9*10E7 |
| Additional remarks: | Storage at 2 - 8°C for 6 months |