

# micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4  
D-18119 Rostock  
Germany

fon: +49 381 / 54 34 56 10  
fax: +49 381 / 54 34 56 20  
e-mail: info@micromod.de



## Technical Datasheet

---

Product-No.:	<b>01-48-103</b>
Product-Name:	<b>micromer ®</b>
surface:	<b>Ni-NTA</b>
size:	<b>1 µm</b>
solid content:	<b>50 mg/ml</b>
quantity:	<b>10 ml</b>
Polydispersity index:	< 0.2
Shape:	spherical
Density:	1.03 g/ccm
Stable in:	aqueous buffers, methanol, ethanol, DMSO
Not stable in:	halogenated hydrocarbons, toluene, strong acidic solutions, eg 10% HCl
Product form:	aqueous suspension
Particles per ml:	9.3*10E10
Particles per mg:	1.9*10E9
Additional remarks:	Storage at 2 - 8°C.