

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

Product-No.:	09-19-252
Product-Name:	nanomag® -D
surface:	streptavidin
size:	250 nm
solid content:	10 mg/ml
quantity:	1 ml
Polydispersity index:	< 0.2
Shape:	cluster-typed
Density:	2.5 g/ccm
Magnetization:	43 emu/g particles (H =1000 Oe)
Saturation magnetization:	>67 emu/g particles (H > 10.000 Oe)
Stable in:	aqueous buffers pH > 4
Not stable in:	organic solvents, acidic solutions pH < 4
Product form:	suspension in PBS buffer (pH = 7.4), stabilized in 0.02 % sodium azide
Particles per ml:	4.9*10E11
Particles per mg:	4.9*10E10
Protein binding capacity:	1.5 - 2.0 µg streptavidin / mg particle (450 - 500 molecules streptavidin per particle)
Additional remarks:	Storage at 2 - 8°C for 3 months