

# micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4  
D-18119 Rostock  
Germany

fon: +49 381 / 54 34 56 10  
fax: +49 381 / 54 34 56 20  
e-mail: info@micromod.de



## Technical Datasheet

---

Product-No.:	<b>29-02-251</b>
Product-Name:	<b>micromer® -greenF</b>
surface:	<b>COOH</b>
size:	<b>25 nm</b>
solid content:	<b>10 mg/ml</b>
quantity:	<b>5 ml</b>
Polydispersity index:	< 0.2
Shape:	spherical
Density:	1.03 g/ccm
Stable in:	aqueous buffers, methanol, ethanol, DMSO
Not stable in:	halogenated hydrocarbons, toluene, strong acidic solutions, eg 10% HCl
Product form:	aqueous suspension
Particles per ml:	1.2*10E15
Particles per mg:	1.2*10E14
Protein binding capacity:	300 - 800 ng protein / mg particle
Optical properties:	green fluorescent
Colour:	yellow
Excitation:	475 nm
Emission:	510 nm
Additional remarks:	Storage at 2 - 8°C.