

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

Product-No.:	30-02-252
Product-Name:	micromer® -redF
surface:	COOH
size:	250 nm
solid content:	25 mg/ml
quantity:	5 ml
Polydispersity index:	< 0.2
Shape:	spherical
Density:	1.03 g/ccm
Stable in:	aqueous buffers, methanol, ethanol, DMSO
Not stable in:	halogenated hydrocarbons, toluene, strong acidic solutions, eg 10% HCl
Product form:	aqueous suspension
Particles per ml:	3.0*10E12
Particles per mg:	1.2*10E11
Protein binding capacity:	300 - 800 ng protein / mg particle
Optical properties:	red fluorescent
Colour:	pink
Excitation:	552 nm
Emission:	580 nm
Additional remarks:	Storage at 2 - 8°C.