

# micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4  
D-18119 Rostock  
Germany

fon: +49 381 / 54 34 56 10  
fax: +49 381 / 54 34 56 20  
e-mail: info@micromod.de



## Technical Datasheet

---

Product-No.:	<b>30-02-603</b>
Product-Name:	<b>micromer® -redF</b>
surface:	<b>COOH</b>
size:	<b>6 µm</b>
solid content:	<b>25 mg/ml</b>
quantity:	<b>10 ml</b>
Standard deviation:	< 5 % (C.V.)
Shape:	spherical
Density:	1.03 g/ccm
Stable in:	aqueous buffers, methanol, ethanol, DMSO
Not stable in:	halogenated hydrocarbons, toluene, strong acidic solutions, eg 10% HCl
Product form:	aqueous suspension
Particles per ml:	2.2*10E8
Particles per mg:	8.6*10E6
Optical properties:	red fluorescent
Colour:	pink
Excitation:	552 nm
Emission:	580 nm
Additional remarks:	Storage at room temperature.