

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

Product-No.:	43-02-204
Product-Name:	sicastar®
surface:	COOH
size:	20 µm
solid content:	50 mg/ml
quantity:	10 ml
Composition:	amorphous silica particles
Standard deviation:	< 11.5 µm
Shape:	spherical with irregularities
Density:	1.8 g/ccm
Porosity:	porous
Stable in:	aqueous buffers, organic solvents
Not stable in:	hydrofluoric acid, strong basic media, e.g. 6 M NaOH
Product form:	aqueous suspension
Particles per ml:	6.6*10E6
Particles per mg:	1.3*10E5
Surface charge density:	4 µmol/g
Protein binding capacity:	540 ng albumin (BSA) / mg particle (357 ng BSA / cm ²)
Additional remarks:	Storage at room temperature.