



micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

Phone: +49 381 / 54 34 56 10
Fax: +49 381 / 54 34 56 20
E-mail Address: info@micromod.de

Technical Data Sheet

Product-No.:	51-00-502
Product-Name:	PLA-Particles -greenF
Surface:	plain
Size:	500 nm
Solid content:	10 mg/ml
Quantity:	5 ml
Composition:	green fluorescent poly(lactic acid) particles
Polydispersity Index:	< 0.2
Shape:	spherical
Density:	1.0 g/ccm
Stable in:	aqueous buffers pH > 4
Not stable in:	organic solvents, acidic solutions pH < 4
Product Form:	aqueous suspension
Particles per ml:	1.5*10E11
Particles per mg:	1.5*10E10
Optical Properties:	green fluorescent
Colour:	yellow
Excitation:	475 nm
Emission:	510 nm
Additional Remarks:	Storage at 2 - 8°C for 6 months.