

# Product Datasheet

## HMGCL Antibody NBP1-32767

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-32767](http://www.novusbio.com/NBP1-32767)

Updated 3/5/2018 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-32767](http://www.novusbio.com/reviews/destination/NBP1-32767)



**NBP1-32767**

HMGCL Antibody


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	0.1M Tris (pH 7.0), 0.1M Glycine and 10% Glycerol
Target Molecular Weight	34 kDa

Product Description	
Host	Rabbit
Gene ID	3155
Gene Symbol	HMGCL
Species	Human, Mouse, Rat
Reactivity Notes	Bovine (88%).
Immunogen	Full length human HMGCL Recombinant protein.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

**Images**


Western Blot: HMGCL Antibody [NBP1-32767] - rat liver lysate/extract 10 % SDS-PAGE gel, antibody dilution 1:10000.




Immunocytochemistry/Immunofluorescence: HMGCL Antibody [NBP1-32767] - HeLa cells were fixed in ice-cold MeOH for 5 min. Green: HMGCL protein stained by HMGCL antibody diluted at 1:500. Blue: Hoechst 33342 staining. Scale bar = 10  $\mu$ m.




Immunohistochemistry-Paraffin: HMGCL Antibody [NBP1-32767] - Paraffin-embedded SW480 xenograft, using antibody at 1:500 dilution.



Western Blot: HMGCL Antibody [NBP1-32767] - (30  $\mu$ g of whole cell lysate) A: H1299 10% SDS PAGE; antibody diluted at 1:3000



Western Blot: HMGCL Antibody [NBP1-32767] - mouse liver lysate/extract 10 % SDS-PAGE gel, antibody dilution 1:10000.





### Novus Biologicals USA

8100 Southpark Way, A-8  
Littleton, CO 80120  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### Novus Biologicals Canada

461 North Service Road West, Unit B37  
Oakville, ON L6M 2V5  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada@novusbio.com

### Novus Biologicals Europe

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info@bio-techne.com

### General Contact Information

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### Products Related to NBP1-32767

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB7156	Goat anti-Rabbit IgG (H+L) Secondary Antibody
NBP2-24891	Rabbit IgG Isotype Control
NBP1-83251PEP	HMGCL Recombinant Protein Antigen

---

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-32767](http://www.novusbio.com/reviews/submit/NBP1-32767)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

