

Product Data Sheet



Product number **A113**
Revision number **RN1.0**

Product Name	(ATTO-550™)cadaverine
Quantity	1 mg
Application	Fluorescent amine donor substrate for transglutaminases
Spectral properties	$\lambda_{Ex, max} = 554 \text{ nm}$, $\lambda_{Em, max} = 576 \text{ nm}$
Molecular Weight	778.4 (perchlorate)
Chemical Structure	Not disclosed
Purity by HPLC	>95 % (254; 550 nm)
Solubility	to be determined
Appearance	Purple solid
Storage	Store at -20°C, desiccate, protect from light
Related products	T001 Bacterial transglutaminase highly purified C076 Z-Gln-Gly-CAD-(ATTO-550™) L106 Atto-550™ TGase Protein Q-Labeling Kit
Reference(s)	None yet
Release date	13 May 2015
NOTE	INTENDED FOR RESEARCH USE ONLY, NOT FOR USE IN HUMAN, THERAPEUTIC OR DIAGNOSTIC APPLICATIONS.