



OXFORD BIOMEDICAL RESEARCH  
P.O. Box 522, Oxford MI 48371 • USA  
USA: 800-692-4633 • Fax: 248-852-4466  
www.oxfordbiomed.com

## Product Specification

### **Anti-Human uPA (FITC labeled)**

**Product #: PA 68**

**Lot #: 400**

Physical Appearance: yellow frozen solution

Solubility: > 10 mg/mL

IgG concentration: 3.2 mg/mL

Labeling Ratio: ~1:1

#### Applications:

Western Blotting: 1:1000 dilution

Elisa: not determined

IHC: not determined

#### Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr

**AVOID REPEATED FREEZE THAW CYCLES**

Buffer: 0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6

Other: < 98% pure by SDS PAGE