



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Rabbit Anti-Mouse PAI-1 (IgG fraction/FITC labeled)

Product #: PI 38

Lot #: ASMPAI-GF-FITC

Molecular Weight: 160 KDa

Physical Appearance: colorless frozen solution

Solubility: >10 mg/mL

Applications:

ELISA: 1 µg/mL

Western Blot: 10 µg/mL

Ultraviolet: Absorbance (280nm) = 5.0 $^{0.1\%} = 1.36$

Concentration: 3.7 mg/mL

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr

AVOID REPEATED FREEZE THAW CYCLES

Buffer: 0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 7.4