



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse tPA Fluorescein Labeled IgG

Product Number: PT 33

Aliquot: 1.0 mg

Lot Number pt33.1105A

Storage: -80°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is fluorescein labeled.
MOLECULAR WEIGHT:	70 kDa
CONCENTRATION:	6.5 mg/mL determined spectrophotometrically.
LABELING RATIO:	0.97 mole FITC/ mole IgG
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS. Protect from light.
IMMUNOGEN:	Recombinant Mouse tPA
CLASS:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
APPLICATIONS:	Western: 5 µg/mL
SPECIES SPECIFICITY:	Mouse.