

Product datasheet for **GE100056**

pCRISPRa-Enhancer Vector, Synergistic Transactivation Domains for CRISPRa

Product data:

Product Type: CRISPR Vectors


Function: CRISPRa Enhancer Vector

Features:

- CRISPRa enhancer vector to robust activate gene expression
- CMV driven MS2-p65-HSF1, transactivation domains for CRISPRa
- Requires pCas-Guide-CRISPRa (GE100055) for gene activation
- MS2-p65-HSF1 synergistically activate gene expression with dCas9-VP64

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)