

Name: Clone OTI4C5, Anti-DDK (FLAG) Magnetic beads
Product Data Sheet - ANTITAG**Catalog: TA150036****Components:** • Clone OTI4C5, Anti-DDK (FLAG) Magnetic beads (TA150036)**Amount:** 5ml**Immunogen:****Host:****Isotype:****Species Reactivity:****Guaranteed Applications:****Suggested Dilutions:** IP, Protein Purification**Concentration:** Lot dependent; please refer to CoA along with shipment**Buffer:****Purification:****Storage Condition:** Shipped at -20C or with ice packs. Upon delivery store at -20C. Dilute in PBS (pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

Target

Target Name:**Alternative Name:** ANTI DDK, ANTI-DDK, ANTIDDK, FLAG TAG, ANTI FLAG, ANTI-FLAG, ANTI FLAG ANTIBODY, FLAG ANTIBODY, M2 ANTIBODY, FLAG M2, ANTIBODY FOR FLAG, 4C5, ANTI DDK 4C5, DDK 4C5, ANTI-DDK 4C5, DYKDDDDK TAG, D-TAG ANTIBODY, ANTI D-TAG, ANTI D TAG, ANTI-D TAG, ANTI-MDDK, ANTI-MYC DDK, ANTI-MYC/DDK, MYC/DDK ANTIBODY, MYC DDK ANTIBODY, ANTI MYC DDK, ANTI DDK MYC, ANTI-DDK/MYC, ANTI-DDK MYC, FLAG AB, FLAG MAB, DDK ANTIBODY, FLAG MAGNETIC, ANTI FLAG MAGNETIC, ANTI FLAG BEAD, ANTI FLAG RESIN, ANTI-FLAG BEAD, ANTI-FLAG RESIN, ANTI-DDK RESIN, ANTI-DDK BEAD, ANTI-DDK COLUMN, ANTI-FLAG COLUMN**Database Link:** -**Function:**

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.

Validation Data

* More validation images may be available on our website: <http://www.origene.com/antibody/TA150036.aspx>

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.

© 2015 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850
