

Name: Recombinant protein of human ficolin (collagen/fibrinogen domain containing) 1 (FCN1)

Catalog: TP304637

Product Data Sheet

Gene Name: Homo sapiens ficolin 1 (FCN1)
GenBank accession: NM_002003

Description: Recombinant protein of human ficolin (collagen/fibrinogen domain containing) 1 (FCN1)

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC204637](#), encoding the full-length human FCN1 with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_001994

Gene Synonym: FCNM; FCNM


Predicted molecular weight: 32.2 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method

Protein Gel Analysis Data:



Purified recombinant protein FCN1 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.