

Name: Recombinant protein of human albumin (ALB)

Catalog: TP305943

Product Data Sheet

Gene Name: Homo sapiens albumin (ALB)
GenBank accession: NM_000477

Description: Recombinant protein of human albumin (ALB)

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC205943](#), encoding the full-length human ALB with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_000468

Gene Synonym: HSA; PRO0883; PRO0903; PRO1341


Predicted molecular weight: 67.2 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method

Protein Gel Analysis Data:



Purified recombinant protein ALB was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

This product is for research and development use only. Not for use in human.