

**Name:** Recombinant protein of human axin 1 (AXIN1), transcript variant 2

**Catalog:** TP308349

## Product Data Sheet

Gene Name: Homo sapiens axin 1 (AXIN1), transcript variant 2  
GenBank accession: NM\_181050

**Description:** Recombinant protein of human axin 1 (AXIN1), transcript variant 2

### Protein Gel Analysis Data:

**Amount:** 20 ug

**Buffer and Storage:** 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

**Source:** Recombinant protein was produced with TrueORF clone, [RC208349](#), encoding the full-length human AXIN1 with C-terminal DDK tag, from human HEK293 cells.

**Protein Accn:** NP\_851393


**Gene Synonym:** AXIN; PPP1R49

**Predicted molecular weight:** 91.5 kDa

**Tags:** C-terminal MYC/DDK

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining.

**Concentration:** >50 ug/mL as determined by microplate BCA method



Purified recombinant protein AXIN1 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

***This product is for research and development use only. Not for use in human.***