

Name: Recombinant protein of human coagulation factor II (thrombin) (F2)

Catalog: TP308589

Product Data Sheet

Gene Name: Homo sapiens coagulation factor II, thrombin (F2), transcript variant 1
GenBank accession: NM_000506

Description: Recombinant protein of human coagulation factor II (thrombin) (F2)

Protein Gel Analysis Data:

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC208589](#), encoding the full-length human F2 with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_000497


Gene Synonym: PT; RPRGL2; THPH1

Predicted molecular weight: 67.5 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method



Purified recombinant protein F2 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

This product is for research and development use only. Not for use in human.

© 2017 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850