

Name: Recombinant protein of human catenin (cadherin-associated protein), beta 1, 88kDa (CTNNB1)

Catalog: TP308947

Product Data Sheet

Gene Name: Homo sapiens catenin beta 1 (CTNNB1), transcript variant 1
GenBank accession: NM_001904

Description: Recombinant protein of human catenin (cadherin-associated protein), beta 1, 88kDa (CTNNB1)

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC208947](#), encoding the full-length human CTNNB1 with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_001895

Gene Synonym: armadillo; CTNNB; MRD19


Predicted molecular weight: 85.3 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method

Protein Gel Analysis Data:



Purified recombinant protein CTNNB1 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.