

**Name:** Recombinant protein of human patatin-like phospholipase domain containing 3 (PNPLA3)

**Catalog:** TP309577

## Product Data Sheet

Gene Name: Homo sapiens patatin like phospholipase domain containing 3 (PNPLA3)  
GenBank accession: NM\_025225

**Description:** Recombinant protein of human patatin-like phospholipase domain containing 3 (PNPLA3)

### Protein Gel Analysis Data:

**Amount:** 20 ug

**Buffer and Storage:** 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

**Source:** Recombinant protein was produced with TrueORF clone, [RC209577](#), encoding the full-length human PNPLA3 with C-terminal DDK tag, from human HEK293 cells.

**Protein Accn:** NP\_079501


**Gene Synonym:** ADPN; C22orf20; iPLA(2)epsilon; iPLA(2)epsilon

**Predicted molecular weight:** 52.7 kDa

**Tags:** C-terminal MYC/DDK

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining.

**Concentration:** >50 ug/mL as determined by microplate BCA method



Purified recombinant protein PNPLA3 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

***This product is for research and development use only. Not for use in human.***

© 2017 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850