

Name: Recombinant protein of human grainyhead-like 1 (Drosophila) (GRHL1), transcript variant 1

Catalog: TP309683

Product Data Sheet

Gene Name: Homo sapiens grainyhead-like 1 (Drosophila) (GRHL1), transcript variant 1
GenBank accession: NM_014552

Description: Recombinant protein of human grainyhead-like 1 (Drosophila) (GRHL1), transcript variant 1

Protein Gel Analysis Data:

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC209683](#), encoding the full-length human GRHL1 with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_055367


Gene Synonym: LBP-32; LBP32; MGR; TF2L2; LBP32; MGR; TF2L2; LBP protein 32; leader-binding protein 32; mammalian grainyhead; transcription factor CP2-like 2; grainyhead-like 1 (Drosophila)

Predicted molecular weight: 69.9 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method



Purified recombinant protein GRHL1 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

This product is for research and development use only. Not for use in human.