

Name: Recombinant protein of human TAR DNA binding protein (TARDBP)

Catalog: TP310639

Product Data Sheet

Gene Name: Homo sapiens TAR DNA binding protein (TARDBP)
GenBank accession: NM_007375

Description: Recombinant protein of human TAR DNA binding protein (TARDBP)

Protein Gel Analysis Data:

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC210639](#), encoding the full-length human TARDBP with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_031401

Gene Synonym: ALS10; TDP-43

Predicted molecular weight: 44.6 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method



Purified recombinant protein TARDBP was analyzed by SDS-PAGE gel and Coomassie Blue Staining.