

Name: Recombinant protein of human peroxisome proliferator-activated receptor gamma (PPARG), transcript variant 3

Catalog: TP312181

Product Data Sheet

Gene Name: Homo sapiens peroxisome proliferator activated receptor gamma (PPARG), transcript variant 3
GenBank accession: NM_138711

Description: Recombinant protein of human peroxisome proliferator-activated receptor gamma (PPARG), transcript variant 3

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC212181](#), encoding the full-length human PPARG with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_619725

Gene Synonym: CIMT1; GLM1; NR1C3; PPARG1; PPARG2; PPARgamma

Predicted molecular weight: 54.5 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method

Protein Gel Analysis Data:



Purified recombinant protein PPARG was analyzed by SDS-PAGE gel and Coomassie Blue Staining.