

Name: Recombinant protein of human autoimmune regulator (AIRE), transcript variant AIRE-1

Catalog: TP313497

Product Data Sheet

Gene Name: Homo sapiens autoimmune regulator (AIRE)
GenBank accession: NM_000383

Description: Recombinant protein of human autoimmune regulator (AIRE), transcript variant AIRE-1

Protein Gel Analysis Data:

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC213497](#), encoding the full-length human AIRE with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_000374


Gene Synonym: AIRE1; APECED; APS1; APSI; PGA1

Predicted molecular weight: 57.5 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method



Purified recombinant protein AIRE was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

This product is for research and development use only. Not for use in human.

© 2017 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850