

**Name:** Purified recombinant protein of Homo sapiens tumor protein p53 (TP53), transcript variant 3

**Catalog:** TP325502

## Product Data Sheet

Gene Name: Homo sapiens tumor protein p53 (TP53), transcript variant 3  
GenBank accession: NM\_001126114

**Description:** Purified recombinant protein of Homo sapiens tumor protein p53 (TP53), transcript variant 3

**Amount:** 20 ug

**Buffer and Storage:** 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

**Source:** Recombinant protein was produced with TrueORF clone, [RC225502](#), encoding the full-length human TP53 with C-terminal DDK tag, from human HEK293 cells.

**Protein Accn:** NP\_001119586

**Gene Synonym:** BCC7; LFS1; P53; TRP53


**Predicted molecular weight:** 37.6 kDa

**Tags:** C-terminal MYC/DDK

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining.

**Concentration:** >50 ug/mL as determined by microplate BCA method

### Protein Gel Analysis Data:



Purified recombinant protein TP53 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

***This product is for research and development use only. Not for use in human.***