

## Product datasheet for **TP500412**

### **Hsd17b14 (NM\_025330) Mouse Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse hydroxysteroid (17-beta) dehydrogenase 14 (Hsd17b14), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Tag:	C-MYC/DDK
Predicted MW:	11.4 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_079606</a></u>
Locus ID:	66065
RefSeq Size:	1365
Cytogenetics:	7 B3
RefSeq ORF:	333
Synonyms:	0610039E24Rik; Dhrrs10; retSDR3

[View online »](#)