

## Product datasheet for **TP503401**

### **Hsd17b10 (NM\_016763) Mouse Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse hydroxysteroid (17-beta) dehydrogenase 10 (Hsd17b10), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Tag:	C-MYC/DDK
Predicted MW:	27.3 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_058043</a></u>
Locus ID:	15108
RefSeq Size:	938
Cytogenetics:	X F3
RefSeq ORF:	786
Synonyms:	17bHSD10; Ads9; ERAB; Hadh2

[View online »](#)