

Name: Purified recombinant protein of Human cytotoxic T-lymphocyte-associated protein 4 (CTLA4), transcript variant 1, with C-terminal DDK/His tag, expressed in human cells Product Data Sheet - PROTEIN	Catalog: TP700216
--	--------------------------

Gene Name: cytotoxic T-lymphocyte associated protein 4
GeneBank accession: NM_005214

Description: Purified recombinant protein of Human cytotoxic T-lymphocyte-associated protein 4 (CTLA4), transcript variant 1, with C-terminal DDK/His tag, expressed in human cells (TP700216)

Amount: 20 ug

Buffer and Storage: 1 x PBS, pH7.4, 10% glycerol

Source: Recombinant protein was produced with A DNA sequence from TrueORF clone, RC210150, encoding the extracellular domain (Ala37-Phe162) of human CTLA4

Protein Accn: NP_005205

Gene Synonym: ALPS5 | CD | CD152 | CELIAC3 | CTLA-4 | GRD4 | GSE | IDDM12


Predicted molecular weight: 16 kDa

Tags: C-terminal DDK/HIS

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Concentration: >50 ug/mL as determined by microplate BCA method

Protein Gel Analysis Data:



Purified recombinant protein CTLA4 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

* More data may be available on our website: [TP700216](#)