

<b>Name: Recombinant protein of human lipopolysaccharide binding protein (LBP)</b> <b>Product Data Sheet - PROTEIN</b>	<b>Catalog: TP720319</b>
---	--------------------------

**Gene Name:** lipopolysaccharide binding protein

**GeneBank accession:** NM\_004139

**Description:** Recombinant protein of human lipopolysaccharide binding protein (LBP) (TP720319)

**Amount:** 10 ug

**Buffer and Storage:** Lyophilized from a 0.2 μM filtered solution of 20mM Phosphate buffer, 150mM NaCl, 1mM EDTA, pH 7.2. Stable for at least 6 months from date of receipt under proper storage and handling conditions.

**Source:** Recombinant protein was produced with A DNA sequence from TrueORF clone, RC221961, encoding the secreted form (Ala26-Val481) of human LBP

**Protein Accn:** NP\_004130

**Gene Synonym:** BPIFD2


**Predicted molecular weight:** 51.9 kDa

**Tags:** C-terminal HIS

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

**Concentration:** Resuspend the protein to the desired concentration in proper buffer.

**Protein Gel Analysis Data:**



Purified recombinant protein LBP was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

\* More data may be available on our website: [TP720319](#)