

Name: Purified recombinant protein of Mouse tumor necrosis factor (Tnf) Product Data Sheet - PROTEIN	Catalog: TP723743
---	--------------------------

Gene Name:	tumor necrosis factor
GeneBank accession:	NM_013693
<hr/>	
Description:	Purified recombinant protein of Mouse tumor necrosis factor (Tnf) (TP723743)
Amount:	10 ug
Buffer and Storage:	10mM NaH ₂ PO ₄ , 150mM NaCl, pH7.2
Source:	Recombinant protein was produced with Mouse TNF-alpha;, the region of Leu80-Leu235, from gene Accession# NM_013693
Protein Accn:	NP_038721
Gene Synonym:	DIF TNF-a TNF-alpha TNFSF2 TNFalpha Tnfa Tnfsf1a Tnlg1f
Predicted molecular weight:	17.3 kDa
Tags:	Tag-free
Purity:	Purity is >98%, as determined by Coomassie stained SDS-PAGE.
Concentration:	Protein formulated and vialled at >50 ug/mL in specified buffer

* More data may be available on our website: [TP723743](#)