

Name: Purified recombinant protein of Mouse interleukin 6 (IL6)

Catalog: TP723747

Product Data Sheet - PROTEIN

Gene Name: interleukin 6

GeneBank accession: NM_031168

Description: Purified recombinant protein of Mouse interleukin 6 (IL6) (TP723747)

Amount: 10 ug

Buffer and Storage: Sodium Acetate and EDTA

Source: Recombinant protein was produced with Mouse IL-6, the region of Phe25-Thr211, from gene Accession# NM_031168

Protein Accn: NP_112445

Gene Synonym: Il-6

Predicted molecular weight: 21.9 kDa

Tags: Tag-free

Purity: Purity is >98%, as determined by Coomassie stained SDS-PAGE.

Concentration: Protein formulated and vialled at >50 ug/mL in specified buffer

* More data may be available on our website: [TP723747](#)