

Name: Purified recombinant protein of Mouse cytotoxic T-lymphocyte-associated protein 4 (Ctla4) Product Data Sheet - PROTEIN	Catalog: TP723873
---	--------------------------

Gene Name: cytotoxic T-lymphocyte-associated protein 4
GeneBank accession: NM_009843

Description: Purified recombinant protein of Mouse cytotoxic T-lymphocyte-associated protein 4 (Ctla4) (TP723873)

Amount: 10 ug

Buffer and Storage: 1 x PBS

Source: Recombinant protein was produced with Mouse CTLA-4, the region of Ala37-Phe162, from gene Accession# NM_009843

Protein Accn: NP_033973

Gene Synonym: Cd152|Ctla-4|Ly-56


Predicted molecular weight: 41 kDa

Tags: Fc chimera

Purity: >95%, as determined by Coomassie stained SDS-PAGE

Concentration: Protein formulated and vialled at >50 ug/mL in specified buffer

Protein Gel Analysis Data:



Purified recombinant protein Ctla4 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

* More data may be available on our website: [TP723873](https://www.origene.com/TP723873)